



Approved for use on (SOMEC 11/03)

DACT-Reg-Me-DMC-R-002-04

Date: 04.12.2003

Report No.: AB-073/03

Page: 01

Department: Life Science & Electronicschemicals  
Analytical Services, Site Elfrut  
Author: DI S. Bär

## REPORT OF ANALYSIS

### Investigation:

Determination of the assay of 4-Pentoxyl-4-terphenylcarboxylic acid in the sample WV 03-157-1 produced in accordance with part C of the Balkevec reference USP 5,948,753 by liquid chromatography (HPLC)

### Sample:

WV 03-157-1

### Customer:

Dr. Nerdinger

### Method:

LM-LC-278

### Raw data:

LC-894/03

### Analytical result:

The current investigation should determine the assay of 4-Pentoxyl-4-terphenylcarboxylic acid in the sample WV 03-157-1 produced in accordance with part C of the Balkevec reference USP 5,948,753.

For the analysis of 4-Pentoxyl-4-terphenylcarboxylic acid a liquid chromatographical method is used. The analysis is carried out by HPLC on a reversed phase silica gel as the stationary phase (elution: acetonitrile/water). UV absorption at three wavelengths (240nm, 270nm and 300nm) is used for detection to confirm the obtained results. The quantification is made by reference material in weight-percent.

The HPLC-analysis leads to nearly the same results at three wavelengths, see table below.

wavelength	Assay of 4-Pentoxyl-4-terphenylcarboxylic acid
240 nm	0.29 % [w/w]
270 nm	0.33 % [w/w]
300 nm	0.25 % [w/w]

The documentation of the HPLC-method is deposited in the raw data and here detailed described. In the chromatograms (Appendix 1) the peaks of interest are labeled with the product name 4-Pentoxyl-4-terphenylcarboxylic acid.

(Stephanie Bär)

Appendix 1

Dr. S. Nerdinger



Date: 04.12.2003  
Report No.: AB-073/03  
Appendix 1

Chromatograms of the sample WV 03-157-1 at three different wavelengths:

